



DMAMUX internal peripheral



Contents

1. DMAMUX internal peripheral	3
-------------------------------------	---



Search for revisions From year (and earlier): From month (and earlier): Tag filter:

Diff selection: Mark the radio boxes of the revisions to compare and hit enter or the button at the bottom.

Legend: **(cur)** = difference with latest revision, **(prev)** = difference with preceding revision, **m** = minor edit.

- **(cur | prev)15:34, 4 February 2020 . . (2,979 bytes) (-118) [Reviewed: 15:42, 4 February 2020 | | [Expert: Aproved, Technical writer: Approved, Maintainer: Approved]]**
- **(cur | prev)16:04, 18 January 2019 m . . (3,097 bytes) (+3,097) [Reviewed: 11:21, 11 February 2019 | | [Expert: Aproved, Technical writer: Approved, Maintainer: Approved]]**

Stable: 04.02.2020 - 15:42 / Revision: 04.02.2020 - 15:34

Search for revisions From year (and earlier): From month (and earlier): Tag filter:

Diff selection: Mark the radio boxes of the revisions to compare and hit enter or the button at the bottom.

Legend: **(cur)** = difference with latest revision, **(prev)** = difference with preceding revision, **m** = minor edit.

- **(cur | prev)15:34, 4 February 2020 . . (2,979 bytes) (-118) [Reviewed: 15:42, 4 February 2020 | | [Expert: Aproved, Technical writer: Approved, Maintainer: Approved]]**
- **(cur | prev)16:04, 18 January 2019 m . . (3,097 bytes) (+3,097) [Reviewed: 11:21, 11 February 2019 | | [Expert: Aproved, Technical writer: Approved, Maintainer: Approved]]**